

Human MUSK (aa 433-783, His & GST Tag) recombinant protein



Catalog Number: 502004

General Information

Gene Name Synonym

Muscle-specific tyrosine-protein kinase receptor

Protein Construction

A DNA sequence encoding the C-terminal segment of human MUSK isoform 2 (O15146-2) (Arg 433-Val 783) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Kinase activity untested

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human MUSK (aa 433-783)/GST chimera consists of 588 amino acids and has a calculated molecular mass of 68 kDa. It migrates as an approximately 58 kDa band in SDS-PAGE under reducing conditions.

Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, pH 7.4, 10mM GSH

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human MUSK Kinase Protein (aa 433-783, His & GST Tag) SDS-PAGE