

Anti-GPR174 antibody

Catalog Number: 176004

Product name

Anti-GPR174 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to GPR174

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human gpr174. aa range:126-175

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.: Western blot analysis of the lysates from HeLa cells using GPR174 antibody.