

# Claudin 9-specific antibody



Catalog Number: 107004

## Product name

Claudin 9-specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Claudin 9-specific Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 25 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_020982). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109343(CLDN9-specific antibody) at dilution of 1:100



IP Result of anti-CLDN9-specific (IP:Catalog No:109343, 4ug; Detection:Catalog No:109343 1:300) with HeLa cells lysate 880ug.