

ZNF557 antibody



Catalog Number: 117004

Product name

ZNF557 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF557 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HepG2 cells, A549 cells, COLO 320 cells, HeLa cells. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of ZNF557 recombinant protein (Accession Number: BC003020). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117004(ZNF557 Antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells using Catalog No:117004(ZNF557 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)