

CCL28 antibody



Catalog Number: 109004

Product name

CCL28 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CCL28 Rabbit Polyclonal antibody. Positive IHC detected in human colon tissue, human colon cancer tissue. Positive WB detected in human spleen tissue, Jurkat cells, PC-3 cells. Observed molecular weight by Western-blot: 8kd-9kd

Preparation

This antibody was obtained by immunization of CCL28 recombinant protein (Accession Number: NM_001301874). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:109004(CCL28 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human colon using Catalog No:109004(CCL28 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon using Catalog No:109004(CCL28 antibody) at dilution of 1:50 (under 40x lens)