

# ERBB4 antibody

Catalog Number: 110404

## Product name

ERBB4 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ERBB4 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 180 kDa

## Preparation

This antibody was obtained by immunization of ERBB4 recombinant protein (Accession Number: NM\_005235). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:500-1:5000

## Validations



IP Result of anti-ERBB4 (IP:Catalog No:110404, 5ug; Detection:Catalog No:110404 1:1000) with mouse brain tissue lysate 3600ug.