

TRPA1 antibody



Catalog Number: 116404

Product name

TRPA1 antibody

Specificity

human; other species not tested.

Antibody description

TRPA1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 120-130kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_007332). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:116404(TRPA1 antibody) at dilution of 1:500