

# INO80 antibody

Catalog Number: 111804

## Product name

INO80 antibody

## Immunogen

[Human INO80 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

INO80 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 177,150kDa

## Preparation

This antibody was obtained by immunization of INO80 recombinant protein (Accession Number: NM\_017553). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111804(INO80 Antibody) at dilution of 1:600



IP Result of anti-INO80 (IP:Catalog No:111804, 4ug; Detection:Catalog No:111804 1:500) with HeLa cells lysate 2000ug.