

NEU1 Polyclonal Antibody

Catalog Number: 164804

Product name

NEU1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to NEU1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 156-415 of human NEU1 (NP_000425.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - NEU1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NEU1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.



Immunohistochemistry - NEU1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using NEU1 antibody at dilution of 1:100



(40x lens).


NEU1 Polyclonal Antibody



Catalog Number: 164804

Immunohistochemistry - NEU1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse spinal cord using NEU1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - NEU1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using NEU1 antibody at dilution of 1:100 (40x lens).