

RPGRIP1 antibody



Catalog Number: 114804

Product name

RPGRIP1 antibody

Immunogen

[Human RPGRIP1 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

RPGRIP1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 70 kDa, 140-150 kDa

Preparation

This antibody was obtained by immunization of RPGRIP1 recombinant protein (Accession Number: NM_001377523). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114804(RPGRIP1 antibody) at dilution of 1:400