

VSX2 antibody



Catalog Number: 116804

Product name

VSX2 antibody

Immunogen

[Human VSX2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

VSX2 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of VSX2 recombinant protein (Accession Number: NM_182894). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:116804(VSX2 Antibody) at dilution of 1:600