

Anti-Acetyl-Histone H3 (Lys27) Rabbit antibody



Catalog Number: 167804

Product name

Monoclonal

Anti-Acetyl-Histone H3 (Lys27) Rabbit antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat

Applications

Antibody description

WB, IHC-P, ICC/IF, FC, IP, ChIP

Acetyl-Histone H3 (Lys27) Rabbit monoclonal antibody

Validations

Preparation

Antigen: A synthetic acetyl-peptide corresponding to residues surrounding Lys27 of human Histone H3


Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality



Western blot detection of Histone H3 (acetyl-Lys27) in HeLa,A549,HL-60 cell lysates using Histone H3 (acetyl-Lys27) Rabbit mAb(1:1000 diluted).Predicted band size:15kDa.Observed band size:17kDa.