

DNAJC12 antibody



Catalog Number: 110104

Product name

DNAJC12 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DNAJC12 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in A549 cells, HEK-293 cells. Observed molecular weight by Western-blot: 23 kDa

Preparation

This antibody was obtained by immunization of DNAJC12 recombinant protein (Accession Number: NM_021800). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:110104(DNAJC12 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110104(DNAJC12 antibody) at dilution of 1:50 (under 10x lens)