

KLHL32 antibody



Catalog Number: 112104

Product name

KLHL32 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

KLHL32 Rabbit Polyclonal antibody. Positive IF detected in Hela.

Preparation

This antibody was obtained by immunization of KLHL32 recombinant protein (Accession Number: NM_001286252). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of Hela cells, using KLHL32 antibody Catalog No:112104 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).