

# Anti-Leptin (7E3) Mouse antibody



Catalog Number: 168104

## Product name

Anti-Leptin (7E3) Mouse antibody

## Specificity

Human

## Antibody description

Leptin (7E3) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant full length of human Leptin protein expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1

## Applications

ELISA

## Validations



Indirect ELISA assay for anti-Leptin mouse mAb. Antigen coating concentration: 4ug/ml.