

# Anti-Thioredoxin/TRX (3C1) Mouse antibody



Catalog Number: 169104

## Product name

Monoclonal

Anti-Thioredoxin/TRX (3C1) Mouse antibody

## Ig Type

## Specificity

Mouse IgG1

Tag

## Applications

## Antibody description

WB

Thioredoxin/TRX (3C1) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant Trx-tag fusion protein expressed in E.coli.

## Formulation


Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

## Validations



Western blot detection of Trx-tag in various bacterial lysates, which containing Trx-tagged proteins (Lane A and Lane B) using Trx-tag mouse mAb (1:1000 diluted).