

EXOC6B antibody

Catalog Number: 110504

Product name

EXOC6B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

EXOC6B Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 94 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_015189). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, IF


Dilutions

Recommended Dilution:


IP: 1:200-1:2000

IF: 1:20-1:200

Validations



IP Result of anti-EXoc6B (IP:Catalog No:110504, 4ug; Detection:Catalog No:110504 1:500) with HEK-293 cells lysate 3200ug.



Immunofluorescent analysis of HEK-293 cells using Catalog No:110504(EXoc6B Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)