

HNRNPD antibody

Catalog Number: 111504

Product name

HNRNPD antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HNRNPD Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, mouse liver tissue. Positive IHC detected in human cervical cancer tissue, human placenta tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 38kd and 40kd

Validations




Preparation

This antibody was obtained by immunization of HNRNPD recombinant protein (Accession Number: NM_002138). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111504(HNRNPD antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:111504(HNRNPD antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:500-1:5000




Immunohistochemical of paraffin-embedded

HNRNPD antibody



Catalog Number: 111504

human cervical cancer using Catalog No:111504(HNRNPD antibody) at dilution of 1:100 (under 40x lens) 

HNRNPD antibody Catalog No:111504 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).

Immunofluorescent analysis of HepG2 cells, using