

MBD4 antibody



Catalog Number: 112504

Product name

MBD4 antibody

Specificity

human, mouse; other species not tested.

Antibody description

MBD4 Rabbit Polyclonal antibody. Positive WB detected in human skin tissue, A549 cells. Observed molecular weight by Western-blot: 66kd

Preparation

This antibody was obtained by immunization of MBD4 recombinant protein (Accession Number: NM_001276270). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:112504(MBD4 antibody) at dilution of 1:500