

BIN3 antibody

Catalog Number: 108504

Product name

BIN3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

BIN3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, A375 cells, HeLa cells, mouse testis tissue. Observed molecular weight by Western-blot: 25kd

Preparation

This antibody was obtained by immunization of BIN3 recombinant protein (Accession Number: XM_047421995). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:117149(BIN3 antibody) at dilution of 1:500