

# GCK antibody



Catalog Number: 110904

## Product name

Polyclonal

GCK antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

GCK Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 52 kDa

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of GCK recombinant protein (Accession Number: NM\_000162). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110904(GCK Antibody) at dilution of 1:600