

AFAP1L2 antibody

Catalog Number: 107904

Product name

AFAP1L2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AFAP1L2 Rabbit Polyclonal antibody. Positive IP detected in mouse thymus tissue. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of AFAP1L2 recombinant protein (Accession Number: NM_001001936). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:107904(AFAP1L2 Antibody) at dilution of



IP Result of anti-AFAP1L2 (IP:Catalog No:107904, 4ug; Detection:Catalog No:107904 1:500) with mouse thymus tissue lysate 2800ug.