

# EFCAB4B antibody



Catalog Number: 110204

## Product name

EFCAB4B antibody

## Specificity

Human; other species not tested.

## Antibody description

EFCAB4B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 46kd

## Preparation

This antibody was obtained by immunization of EFCAB4B recombinant protein (Accession Number: NM\_032680). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110204(EFCAB4B Antibody) at dilution of 1:600