

GRHPR antibody

Catalog Number: 111204

Product name

GRHPR antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GRHPR Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human liver tissue, L02 cells, MCF-7 cells. Positive IHC detected in human liver tissue, human ovary tumor tissue. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of GRHPR recombinant protein (Accession Number: NM_012203). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111204(GRHPR Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:111204(GRHPR Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:111204(GRHPR Antibody) at dilution of 1:200

GRHPR antibody



Catalog Number: 111204

(under 40x lens)