

LGALS3BP antibody



Catalog Number: 112204

Product name

LGALS3BP antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

LGALS3BP Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse brain tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human breast cancer tissue, human skin tissue. Observed molecular weight by Western-blot: 65-70 kDa




Preparation

This antibody was obtained by immunization of LGALS3BP recombinant protein (Accession Number: NM_005567). Purification method: Antigen Affinity purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112204(LGALS3BP antibody) at dilution of

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



1:1000

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:112204(LGALS3BP Antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:500-1:5000




LGALS3BP antibody



Catalog Number: 112204

Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:112204(LGALS3BP Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:112204(LGALS3BP Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)