

Human CD300A (Fc Tag) recombinant protein



Catalog Number: 503204

General Information

Gene Name Synonym

CD300 antigen-like family member A; CMRF-35-H9; CMRF35-H; IRC1/IRC2; Immunoglobulin superfamily member 12; Inhibitory receptor protein 60; NK inhibitory receptor

Protein Construction

A DNA sequence encoding the human CD300A (Q9UGN4-1)(Met1-Gln178) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

(97.9+1.1) % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 18

Molecular Mass

The recombinant human CD300A/Fc is a disulfide-linked homodimer. The reduced monomer comprises 402 amino acids and has a predicted molecular mass of 44.5 kDa. The apparent molecular mass of the protein is approximately 58-66 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CD300A Protein (Fc Tag) SDS-PAGE