

# NKX2-5 antibody



Catalog Number: 113204

## Product name

NKX2-5 antibody

## Specificity

Human; other species not tested.

## Antibody description

NKX2-5 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 37kd

## Preparation

This antibody was obtained by immunization of NKX2-5 recombinant protein (Accession Number: NM\_004387). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113204(NKX2-5 antibody) at dilution of 1:1200



Immunohistochemical of paraffin-embedded human brain using Catalog No:113204(NKX2-5 antibody) at dilution of 1:100 (under 10x lens)