

TMIE antibody

Catalog Number: 116204

Product name

TMIE antibody

Immunogen

[Human TMIE Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

TMIE Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of TMIE recombinant protein (Accession Number: NM_147196). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116204(TMIE Antibody) at dilution of 1:1000