

# Anti-RhoGDI Rabbit antibody



Catalog Number: 167604

## Product name

Anti-RhoGDI Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

RhoGDI Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human RhoGDI

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions

WB: 1/1000

ICC/IF: 1/20

## Validations



Western blot detection of RhoGDI in  
HeLa,A549,HL-60,C6 using RhoGDI Rabbit  
mAb(1:1000 diluted)