

Human CA11/CARP2/UNQ211/PRO237 (His tag) recombinant protein



Catalog Number: 519604

General Information

Gene Name Synonym

Carbonic Anhydrase-Related Protein 11, CA-RP XI, CA-XI, CARP XI, Carbonic Anhydrase-Related Protein 2, CA-RP II, CARP-2, CA11, CARP2

Protein Construction

A DNA sequence encoding the Human CA11/CARP2/UNQ211/PRO237 His24-Arg328 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human CA11/CARP2/UNQ211/PRO237 (His tag)
recombinant protein