

# HRP-conjugated Mouse anti DDDDK-Tag mAb



Catalog Number: 166304

## Product name

HRP-conjugated Mouse anti DDDDK-Tag mAb

## Antibody description

Monoclonal antibody to HRP-conjugated Mouse anti DDDDK-Tag

## Preparation

Antigen: A synthetic peptide of DDDDK-Tag

## Formulation

PBS, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG


## Applications

WB

## Dilutions

WB 1:2000 - 1:10000

## Validations



Western blot - HRP-conjugated Mouse anti DDDDK-Tag mAb

Western blot analysis of recombinant DDDDK protein with HRP-conjugated Mouse anti DDDDK-Tag mAb at dilution of 1:5000. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.