

# Ch-TOG antibody



Catalog Number: 109304

## Product name

Ch-TOG antibody

## Immunogen

[Human CKAP5 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

Ch-TOG Rabbit Polyclonal antibody. Positive WB detected in NCCIT cell. Observed molecular weight by Western-blot: 220 kDa

## Preparation

This antibody was obtained by immunization of Ch-TOG recombinant protein (Accession Number: NM\_001008938). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



NCCIT cell were subjected to SDS PAGE followed by western blot with Catalog No:109304(CKAP5 Antibody) at dilution of 1:1000