

# VAMP1 antibody

Catalog Number: 116704

## Product name

VAMP1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

VAMP1 Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 13 kDa

## Preparation

This antibody was obtained by immunization of VAMP1 recombinant protein (Accession Number: NM\_199245). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116704(VAMP1 antibody) at dilution of 1:500



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:116704(VAMP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)