

PLSCR4 antibody



Catalog Number: 114014

Product name

PLSCR4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PLSCR4 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human brain tissue, human liver tissue, human placenta tissue. Observed molecular weight by Western-blot: 37kd

Preparation

This antibody was obtained by immunization of PLSCR4 recombinant protein (Accession Number: XM_011513031). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:114014(PLSCR4 antibody) at dilution of 1:400