

Anti-MACO1 antibody

Catalog Number: 175014

Product name

Anti-MACO1 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to MACO1

Preparation

This antigen of this antibody was synthetic peptide within human tmem57 aa 119-171.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:2,000-1:5,000

Validations




Fig1: Western blot analysis on human kidney (A) and D3 (B) using anti- TMEM57 polyclonal antibody.