

SCP2 antibody

Catalog Number: 115014

Product name

SCP2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SCP2 Rabbit Polyclonal antibody. Positive IF detected in BxPC-3 cells. Positive WB detected in A549 cells, BxPC-3 cells, mouse testis tissue, PC-3 cells. Observed molecular weight by Western-blot: 44kd, 52kd

Preparation

This antibody was obtained by immunization of SCP2 recombinant protein (Accession Number: NM_001007098). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:115014(SCPX antibody) at dilution of 1:500



Immunofluorescent analysis of BxPC-3 cells using Catalog No:115014(SCPX Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)