

Tislelizumab reference antibody



Catalog Number: 177014

Product name

Tislelizumab reference antibody

Specificity

Human PDCD1 / PD-1 / CD279

Antibody description

Anti-PDCD1 / PD-1 / CD279 Antibody
(Tislelizumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Glycine, 0.2 NaCl, 50 mM Tris-HCl, pH 6.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG4SP

Applications


ELISA

Validations



SDS-PAGE

Anti-PDCD1 / PD-1 / CD279 Antibody (tislelizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-PDCD1 / PD-1 / CD279 Antibody (tislelizumab) is more than 95%, determined by SEC-HPLC.