

Anti-ACAT1 Rabbit antibody



Catalog Number: 169014

Product name

Anti-ACAT1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ACAT1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human ACAT1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of ACAT1 in K562, HeLa cell lysates using ACAT1 Rabbit mAb(1:1000 diluted). Predicted band size: 45kDa. Observed band size: 45kDa.