

# Human KPNB1/NTF97 (His tag) recombinant protein



Catalog Number: 519014

## General Information

### Gene Name Synonym

Importin subunit beta-1; Importin-90; Karyopherin subunit beta-1; Nuclear factor p97; Pore targeting complex 97 kDa subunit; PTAC97; KPNB1; NTF97

### Protein Construction

A DNA sequence encoding the Human KPNB1/NTF97 Met1-Ala876 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 1mM DTT, 30% glycerol, 0.1M NaCl, pH8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human KPNB1/NTF97 (His tag) recombinant protein