

# RAB17 antibody

Catalog Number: 114414

---

## Product name

RAB17 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RAB17 Rabbit Polyclonal antibody. Positive IF detected in WIF-B cells.

## Preparation

This antibody was obtained by immunization of RAB17 recombinant protein (Accession Number: NM\_022449). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF


## Dilutions

Recommended Dilution:

IF: N/A

## Validations

### rab17 antibody



IF result of RAB17 antibody (Catalog No:111366) with WIF-B cells overexpressed myc-tagged RAB17 by A.C. Striz and Dr. P.L. Tuma.