

SMA4 antibody

Catalog Number: 115414

Product name

SMA4 antibody

Specificity

human; other species not tested.

Antibody description

SMA4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 15kd

Preparation

This antibody was obtained by immunization of SMA4 recombinant protein (Accession Number: BC002622). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115414(SMA4 antibody) at dilution of 1:300