

Phospho-IL10RA-S319/S323 pAb



Catalog Number: 166414

Product name

Phospho-IL10RA-S319/S323 pAb

Specificity

Human

Antibody description

Polyclonal antibody to Phospho-IL10RA-S319/S323

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding S319/S323 of human IL10RA

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-IL10RA-S319/S323 pAb

Western blot analysis of extracts from HeLa cells, using Phospho-IL10RA-S319/S323 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.