

Anti-Cytochrome C Rabbit antibody



Catalog Number: 167414

Product name

Rabbit IgG

Anti-Cytochrome C Rabbit antibody

Applications

Specificity

WB, ICC/IF, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

Cytochrome C Rabbit polyclonal antibody

ICC/IF: 1/20-1/50

Preparation

IP: 1/20

Antigen: A synthetic peptide of human Cytochrome C

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type



Western blot detection of Cytochrome C in K562, Rat Brain, C6, 3T3, HeLa cell lysates using Cytochrome C Rabbit pAb(1:1000 diluted). Predicted band size: 12kDa. Observed band size: 14kDa.