

Anti-CABP1 antibody

Catalog Number: 175814

Product name

Anti-CABP1 antibody

Specificity

Rat, Mouse

Antibody description

Rabbit polyclonal antibody to CABP1

Preparation

This antigen of this antibody was synthetic peptide of cabp

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1:: Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue with CABP Rabbit pAb diluted at 1:2,000.




Fig2:: Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using CABP Rabbit pAb diluted at 1:200.