

STXBP6 antibody

Catalog Number: 115814

Product name

STXBP6 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

STXBP6 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human heart tissue. Positive IHC detected in human lung cancer. Observed molecular weight by Western-blot: 28-30kd

Preparation

This antibody was obtained by immunization of STXBP6 recombinant protein (Accession Number: NM_014178). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115814(STXBP6 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:115814(STXBP6 antibody) at dilution of 1:50 (under 10x lens)