

# Anti-MYO19 Rabbit antibody



Catalog Number: 168814

## Product name

Anti-MYO19 Rabbit antibody

## Specificity

Human

## Antibody description

MYO19 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human MYO19

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions

1:1000

1:20

## Validations



Western blot detection of MYO19 in K562,3T3,HeLa cell lysates using MYO19 Rabbit mAb(1:1000 diluted).Predicted band size:109kDa.Observed band size:109kDa.