

# CADM4 antibody



Catalog Number: 108814

---

## Product name

CADM4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CADM4 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, COLO 320 cells, mouse lung tissue. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of CADM4 recombinant protein (Accession Number: NM\_145296). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108814(CADM4 antibody) at dilution of 1:1000