

Anti-Lamin B2 Rabbit antibody



Catalog Number: 167114

Product name

Anti-Lamin B2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Lamin B2 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Lamin B2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/2000

FC: 1/20

Validations



Western blot detection of Lamin B2 in K562, Rat Brain, C6, 3T3, HeLa cell lysates using Lamin B2 Rabbit pAb(1:1000 diluted). Predicted band size: 68kDa. Observed band size: 68kDa.