

Tinurilimab reference antibody



Catalog Number: 177114

Product name

Tinurilimab reference antibody

Specificity

Human CEACAM6 / CD66c

Antibody description

Anti-CEACAM6 / CD66c Antibody (Tinurilimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG2SA

Applications


ELISA

Validations




SDS-PAGE

Anti-CEACAM6 / CD66c Antibody (tinurilimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-CEACAM6 / CD66c Antibody (tinurilimab) is 95.93%, determined by SEC-HPLC.



Bioactivity: ELISA

Immobilized human CEACAM6 His at 2 µg/mL can bind Anti-CEACAM6 / CD66c Antibody (tinurilimab), EC50=0.005084 µg/mL.