

Mouse CD63 / Tspan-30 / Tetraspanin-30 recombinant protein



Catalog Number: 503514

General Information

Protein Construction

A DNA sequence encoding the rat Cd63 (NP_031679.1) (Ala103-Ile203) was expressed with two additional amino acids (Gly & Pro) at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> (5.6+73.7+8.4) as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly

Molecular Mass

The recombinant rat Cd63 consists 103 amino acids and predicts a molecular mass of 11.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat CD63 / Tspan-30 / Tetraspanin-30 Protein